

DERMACLIP® NON-INVASIVE SKIN CLOSURE, REGULAR CLIPS (62001 / 62002)

Product Description

The DermaClip® Skin Closure Device is a needle-free, non-invasive alternative to sutures, staples and glues. Proven both operationally and clinically, the DermaClip device allows the provider to make a repeatable, quality closure, often in less than a minute, even under adverse conditions. Use of the device is simple and can be learned in minutes.

Compact and lightweight – a box (containing 12 individual pkgs) providing 12" (2-Clip Pkg) to 36" (5-Clip pkg) of closure lasting 7+ days, weighs only 2.3 oz and fits anywhere. The DermaClip device is usable on all types of skin and can hold high tension wounds. DermaClip devices provide the go-to, easy-to-use device for wound closure in the field or in the clinic.

Benefits & Attributes

- **Fast & Strong:** Close wounds in less than a minute
- **Simple & Safe:** Mastered in minutes, no needles, easy removal
- **Multi-Functional:** For temporary or final closure
- **Versatile:** Usable on straight, jagged or curving wounds
- **Stackable:** Use multiple packages for wounds of any length
- **Tactile:** Safely apply in low visibility with tactile feedback
- **Force Multiplier:** Helps keep warfighter on the battlefield
- **72-Hour Sustainability:** Lasts 7+ days
- **All Roles of Care:** Works for wound care in ROC 1, ROC 2, & ROC 3
- **Fits Anywhere:** Compact, only 0.13 oz/pkg; 2.3 oz/12 pkg box

Details


- **Part Number:** 62001 (1 inch/pkg), 62002 (3 inch/pkg)
- **NSN:** 62001 – pending | 62002 – 6515-01-662-6938
- **Unit of Issue U/I:** Box of 12 Packages (1 Pkg = EA)
- **Source of Supply:** Beacon Point Associates (SDVOSB)
UE ID: L97MMG5VW9N7 **CAGE:** 6U5B7
- **DAPA #** SPO200-18-H-0042, **ECAT #** SPE2DE-21-D-7004
- **Service Life:** 3 years
- **Weight U/I:** 62001 – 2.1 oz (EA: 0.10 oz) | 62002 – 2.3 oz (EA: 0.13 oz)
- **Cube: (H x W x D):** Box 4.6" x 5.8" x 1.5" (EA: 4.1" x 5.3" x 0.1")
- **Clip Dimensions (Single Clip):** 0.43" x 2.25"

military@dermaclip.com

Phone: 713-682-3185

DermaClip®

P/N
62001



1 Package, 1" of Closure
1 Box = 12 Pkg with 2 Clips



P/N
62002



1 Package, 3" of Closure
1 Box = 12 Pkg with 5 Clips



DERMACLIP® NON-INVASIVE SKIN CLOSURE, LARGE CLIPS (62003 / 62004)

Product Description

The DermaClip® Skin Closure Device is a needle-free, non-invasive alternative to sutures, staples and glues. Proven both operationally and clinically, the DermaClip device allows the provider to make a repeatable, quality closure, often in less than a minute, even under adverse conditions. Use of the device is simple and can be learned in minutes.

Compact and lightweight – a box (containing 12 individual pkgs) providing 12" to 36" of closure lasting 5+ days, weighs only 2 oz and fits anywhere. DermaClip is usable on all types of skin and can hold high tension wounds. The device provides the go-to, easy-to-use device for wound closure in the field or in the clinic.

Benefits & Attributes

- **Fast & Strong:** Close high-tension wounds in less than a minute
- **Simple & Safe:** Mastered in minutes, no needles, easy removal
- **Multi-Functional:** For temporary or final closure
- **Versatile:** Usable on straight, jagged or curving wounds
- **Stackable:** Use multiple packages for wounds of any length
- **Tactile:** Safely apply in low visibility with tactile feedback
- **Force Multiplier:** Helps keep warfighter on the battlefield
- **72-Hour Sustainability:** Lasts 5+ days
- **All Roles of Care:** Works for wound care in ROC 1, ROC 2, & ROC 3
- **Fits Anywhere:** Compact, only 0.15 oz/pkg; 2.7 oz/12pkg box

Details

- **Part Number:** 62003 (1 inch/pkg), 62004 (3 inch/pkg)
- **NSN:** 62003 – pending | 62004 – 6515-01-661-8487
- **Unit of Issue U/I:** Box of 12 Packages (1 Pkg = EA)
- **Source of Supply:** Beacon Point Associates (SDVOSB)
UE ID: L97MMG5VW9N7 **CAGE:** 6U5B7
- **DAPA #** SPO200-18-H-0042, **ECAT #** SPE2DE-21-D-7004
- **Service Life:** 3 years
- **Weight U/I:** 62003 – 2.0 oz (EA: 0.09 oz) | 62004 – 2.5 oz (EA: 0.14 oz)
- **Cube: (H x W x D):** Box 4.6" x 5.8" x 1.5" (EA: 4.1" x 5.3" x 0.1")
- **Clip Dimensions (Single Clip):** 0.81" x 2.41

military@dermaclip.com

Phone: 713-682-3185

DermaClip®

P/N
62003



1 Package, 1" of Closure
1 Box = 12 Pkg with 1 Clip



P/N
62004



1 Package, 3" of Closure
1 Box = 12 Pkg with 3 Clips



DERMACLIP® EXTREME, NON-INVASIVE SKIN CLOSURE, LARGE CLIPS (62005)

Product Description

The DermaClip® Skin Closure Device is a needle-free, non-invasive alternative to sutures, staples and glues. Proven both operationally and clinically, the DermaClip device allows the provider to make a repeatable, quality closure, often in less than a minute, even under adverse conditions. Use of the device is simple and can be learned in minutes.

Compact and lightweight – a box (containing 12 individual pkgs) providing 12" to 36" of closure lasting 5+ days, weighs only 2 oz and fits anywhere. DermaClip is usable on all types of skin and can hold high tension wounds. The device provides the go-to, easy-to-use device for wound closure in the field or in the clinic.

DermaClip Extreme adhesive is 3.5x stronger than the standard product.

Benefits & Attributes

- **Fast & Strong:** Close high-tension wounds in less than a minute
- **Simple & Safe:** Mastered in minutes, no needles, easy removal
- **Multi-Functional:** For temporary or final closure
- **Versatile:** Usable on straight, jagged or curving wounds
- **Stackable:** Use multiple packages for wounds of any length
- **Tactile:** Safely apply in low visibility with tactile feedback
- **Force Multiplier:** Helps keep warfighter on the battlefield
- **72-Hour Sustainability:** Lasts 5+ days
- **All Roles of Care:** Works for wound care in ROC 1, ROC 2, & ROC 3
- **Fits Anywhere:** Compact, only 0.15 oz/pkg; 2.5 oz/12 pkg box

Details


- **Part Number:** 62005 (3 inch/pkg)
- **NSN:** 6515-01-673-7948
- **Unit of Issue U/I:** Box of 12 Packages (1 Pkg = EA)
- **Source of Supply:** Beacon Point Associates (SDVOSB)
UE ID: L97MMG5VW9N7 **CAGE:** 6U5B7
- **DAPA #** SPO200-18-H-0042, **ECAT #** SPE2DE-21-D-7004
- **Service Life:** 4 years
- **Weight U/I:** 2.7 oz (EA: 0.15 oz)
- **Cube: (H x W x D):** Box 4.6" x 5.8" x 1.5" (EA: 4.1" x 5.3" x 0.1")
- **Clip Dimensions (Single Clip):** 0.81" x 2.41

military@dermaclip.com

Phone: 713-682-3185

DermaClip®

P/N
62005



1 Package, 3" of Closure
1 Box = 12 Pkg with 3 Clips



DERMAQUIK® NON-INVASIVE SKIN CLOSURE (62501 / 62502)

P/N
62501

(Box of 12)



1 Package, 1" of Closure
1 Box = 12 Pkg with 1 Clip



P/N
62502

(Box of 12)



Product Description

The DermaQuik® Skin Closure Device is a needle-free, non-invasive alternative to sutures, staples and glues. The DermaQuik device allows repeatable, quality closure, often in less than a minute, even under adverse conditions. The device can be closed one-handed, allowing for more versatile use in self-care applications.

Compact and lightweight – a box (containing 12 pkgs) providing 12" of closure lasting 5+ days, weighs only 2 oz and fits anywhere. DermaQuik is usable on all types of skin and can hold high tension wounds. The device provides the go-to, easy-to-use device for wound closure in the field or clinic.

The DermaQuik Sport & Extreme adhesive is 3.5x stronger than the Home & Fields product.

Details

- **Part Number:** 62501 (Home & Field), 62502 (Sport & Extreme)
- **NSN:** 62501 – pending | 62502 – pending
- **Unit of Issue U/I:** Box of 12 Packages (1 Pkg = EA)
- **Source of Supply:** Beacon Point Associates (SDVOSB)
UE ID: L97MMG5VW9N7 **CAGE:** 6U5B7
- **DAPA #** SPO200-18-H-0042, **ECAT #** SPE2DE-21-D-7004
- **Service Life:** 62501 - 3 years | 62502 - 4 years
- **Weight U/I:** 62501 – 2.1 oz (EA: 0.1 oz) | 62502 – 2.1 oz (EA: 0.1 oz)
- **Cube: (H x W x D):** Box 4.6" x 5.8" x 1.5" (EA: 4.1" x 5.3" x 0.1")
- **Clip Dimensions (Single Clip):** 0.81" x 2.41

Benefits & Attributes

- **Fast & Strong:** Close high-tension wounds in less than a minute
- **Simple & Safe:** Mastered in minutes, no needles, easy removal
- **Multi-Functional:** For temporary or final closure
- **Versatile:** Usable on straight, jagged or curving wounds; Close with 1 hand
- **Stackable:** Use multiple packages for wounds of any length
- **Tactile:** Safely apply in low visibility with tactile feedback
- **Force Multiplier:** Helps keep warfighter on the battlefield
- **72-Hour Sustainability:** Lasts 5+ days
- **All Roles of Care:** Works for wound care in ROC 1, ROC 2, & ROC 3
- **Fits Anywhere:** Compact, only 0.10 oz/pkg; 2.1 oz/12 pkg box

military@dermaclip.com Phone: 713-682-3185

DermaQuik® by DermaClip®

PRODUCT NUMBER	ITEM	DESCRIPTION	U/M	GOVERNMENT PRICE	NSN
62001	DermaClip Regular Clip, 2 Pack (1"), Standard Adhesive, Box (12 Packages)	DermaClip Regular sized devices (0.43" width per clip) with Standard adhesive, 2 clips per Package (closing 1" of wound per Pkg). Sold in a Box of 12 Packages. DermaClip is a needle-free, non-invasive alternative to sutures, staples and glues that allows the provider to make a repeatable, quality closure, often in less than a minute, even under adverse conditions. Compact and lightweight, DermaClip is usable in any environment, from the field to the surgical suite, on all types of skin and can hold high tension wounds.	BX (12 EA)	\$ 230.00	pending
62002	DermaClip Regular Clip, 5 Pack (3"), Standard Adhesive, Box (12 Packages)	DermaClip Regular sized devices (0.43" width per clip) with Standard adhesive, 5 clips per Package (closing 3" of wound per Pkg). Sold in a Box of 12 Packages. DermaClip is a needle-free, non-invasive alternative to sutures, staples and glues that allows the provider to make a repeatable, quality closure, often in less than a minute, even under adverse conditions. Compact and lightweight, DermaClip is usable in any environment, from the field to the surgical suite, on all types of skin and can hold high tension wounds.	BX (12 EA)	\$ 350.00	6515-01-662-6938
62003	DermaClip Large Clip, 1 Pack (1"), Standard Adhesive, Box (12 Packages)	DermaClip Large sized devices (0.81" width per clip) with Standard adhesive, 1 clip per Package (closing 1" of wound per Pkg). Sold in a Box of 12 Packages. DermaClip is a needle-free, non-invasive alternative to sutures, staples and glues that allows the provider to make a repeatable, quality closure, often in less than a minute, even under adverse conditions. Compact and lightweight, DermaClip is usable in any environment, from the field to the surgical suite, on all types of skin and can hold high tension wounds.	BX (12 EA)	\$ 230.00	pending
62004	DermaClip Large Clip, 3 Pack (3"), Standard Adhesive, Box (12 Packages)	DermaClip Large sized devices (0.81" width per clip) with Standard adhesive, 3 clips per Package (closing 3" of wound per Pkg). Sold in a Box of 12 Packages. DermaClip is a needle-free, non-invasive alternative to sutures, staples and glues that allows the provider to make a repeatable, quality closure, often in less than a minute, even under adverse conditions. Compact and lightweight, DermaClip is usable in any environment, from the field to the surgical suite, on all types of skin and can hold high tension wounds.	BX (12 EA)	\$ 350.00	6515-01-661-8487
62005	DermaClip EXTREME Large Clip, 3 Pack (3"), EXTREME Adhesive, Box (12 Packages)	DermaClip EXTREME Large sized devices (0.81" width per clip) with Extreme adhesive (3.5x strength of Standard adhesive), 3 clips per Package (closing 3" of wound per Pkg). Sold in a Box of 12 Packages. DermaClip EXTREME is a needle-free, non-invasive alternative to sutures, staples and glues that allows the provider to make a repeatable, quality closure, often in less than a minute, even under adverse conditions. Compact and lightweight, DermaClip is usable in any environment, from the field to the surgical suite, on all types of skin and can hold high tension wounds.	BX (12 EA)	\$ 495.00	6515-01-673-7948

PRODUCT NUMBER	ITEM	DESCRIPTION	U/M	GOVERNMENT PRICE	NSN
63001	DermaClip Regular Clip, 2 Pack (1"), Standard Adhesive, Case (12 Boxes, 144 Total Packages)	DermaClip Regular sized devices (0.43" width per clip) with Standard adhesive, 2 clips per Package (closing 1" of wound per Pkg). Sold in a Case of 12 Boxes (each Box contains 12 Packages). DermaClip is a needle-free, non-invasive alternative to sutures, staples and glues that allows the provider to make a repeatable, quality closure, often in less than a minute, even under adverse conditions. Compact and lightweight, DermaClip is usable in any environment, from the field to the surgical suite, on all types of skin and can hold high tension wounds.	CS (12 BXS)	\$ 2,615.00	pending
63002	DermaClip Regular Clip, 5 Pack (3"), Standard Adhesive, Case (12 Boxes, 144 Total Packages)	DermaClip Regular sized devices (0.43" width per clip) with Standard adhesive, 5 clips per Package (closing 3" of wound per Pkg). Sold in a Case of 12 Boxes (each Box contains 12 Packages). DermaClip is a needle-free, non-invasive alternative to sutures, staples and glues that allows the provider to make a repeatable, quality closure, often in less than a minute, even under adverse conditions. Compact and lightweight, DermaClip is usable in any environment, from the field to the surgical suite, on all types of skin and can hold high tension wounds.	BX (12 BXS)	\$ 4,055.00	6515-01-662-6997
63003	DermaClip Large Clip, 1 Pack (1"), Standard Adhesive, Case (12 Boxes, 144 Total Packages)	DermaClip Large sized devices (0.81" width per clip) with Standard adhesive, 1 clip per Package (closing 1" of wound per Pkg). Sold in a Case of 12 Boxes (each Box contains 12 Packages). DermaClip is a needle-free, non-invasive alternative to sutures, staples and glues that allows the provider to make a repeatable, quality closure, often in less than a minute, even under adverse conditions. Compact and lightweight, DermaClip is usable in any environment, from the field to the surgical suite, on all types of skin and can hold high tension wounds.	BX (12 BXS)	\$ 2,615.00	pending
63004	DermaClip Large Clip, 3 Pack (3"), Standard Adhesive, Case (12 Boxes, 144 Total Packages)	DermaClip Large sized devices (0.81" width per clip) with Standard adhesive, 3 clips per Package (closing 3" of wound per Pkg). Sold in a Case of 12 Boxes (each Box contains 12 Packages). DermaClip is a needle-free, non-invasive alternative to sutures, staples and glues that allows the provider to make a repeatable, quality closure, often in less than a minute, even under adverse conditions. Compact and lightweight, DermaClip is usable in any environment, from the field to the surgical suite, on all types of skin and can hold high tension wounds.	CS (12 BXS)	\$4,055.00	6515-01-662-6961
63005	DermaClip EXTREME Large Clip, 3 Pack (3"), EXTREME Adhesive, Case (12 Boxes, 144 Total Packages)	DermaClip EXTREME Large sized devices (0.81" width per clip) with Extreme adhesive (3.5x strength of Standard adhesive), 3 clips per Package (closing 3" of wound per Pkg). Sold in a Case of 12 Boxes (each Box contains 12 Packages). DermaClip EXTREME is a needle-free, non-invasive alternative to sutures, staples and glues that allows the provider to make a repeatable, quality closure, often in less than a minute, even under adverse conditions. Compact and lightweight, DermaClip is usable in any environment, from the field to the surgical suite, on all types of skin and can hold high tension wounds.	CS (12 BXS)	\$ 5,795.00	6515-01-673-7942

PRODUCT NUMBER	ITEM	DESCRIPTION	U/M	GOVERNMENT PRICE	NSN
62501	DermaQuik Home & Field, 1 Pack (1"), Standard Adhesive, Box (12 Packages)	<p>DermaClip Regular sized devices (0.43" width per clip) with Standard adhesive, 2 clips per Package (closing 1" of wound per Pkg). Sold in a Box of 12 Packages.</p> <p>DermaClip is a needle-free, non-invasive alternative to sutures, staples and glues that allows the provider to make a repeatable, quality closure, often in less than a minute, even under adverse conditions. Compact and lightweight, DermaClip is usable in any environment, from the field to the surgical suite, on all types of skin and can hold high tension wounds.</p>	BX (12 EA)	\$ 171.50	pending
62502	DermaQuik Sport & Extreme, 1 Pack (1"), Extreme Adhesive, Box (12 Packages)	<p>DermaClip Regular sized devices (0.43" width per clip) with Standard adhesive, 5 clips per Package (closing 3" of wound per Pkg). Sold in a Box of 12 Packages.</p> <p>DermaClip is a needle-free, non-invasive alternative to sutures, staples and glues that allows the provider to make a repeatable, quality closure, often in less than a minute, even under adverse conditions. Compact and lightweight, DermaClip is usable in any environment, from the field to the surgical suite, on all types of skin and can hold high tension wounds.</p>	BX (12 EA)	\$ 192.50	pending